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规格: 100mL/500mL

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**Cat:** G1100

**Size:** 100mL/500mL

**Storage:** RT, avoid light, valid for 1 year.

## Introduction

Eosin is an acidic dye that interacts with the cytoplasmic proteins rich in basic amino acids to form a complex red-pink brilliant. Eosin Y Stain Solution, Water Solvent is a ready-to-use solution that can be used for routine staining in the hematoxylin-eosin staining process (H&E). The H&E staining is the most frequently used staining method for histological material.

## Protocol(for reference only)

### 1. Sample Treatment

For paraffin section: dewax in xylene twice for 5-10min each. Rehydrate with series ethanol (100%,95%,85%,75%) for each level 3min,and finally in distilled water for 3min.

For frozen section: restore to RT and wash with distilled water for 2min.

For cultured cell: fix with 4% PFA for more than 10min.Then wash with distilled water twice for 2min each..

### 2. Eosin Staining

① Staining with Eosin for 2-5 min (the time can be adjusted according to the staining results and experimental requirements).

② Quickly wash the excess dye with distilled water for 2-3s.

③ Dehydrate with series ethanol(75%,85%,95%) for 2-3s each.

④ 100% ethanol (I) for 2-3s,100% ethanol (II) for 1 min.

⑤ Wash with xylene twice for 1 min each.

⑥ Seal with resinene and view under the microscope.

*Note: if it is used for re-staining after immunohistochemical staining, you can refer to the above steps to directly carry out eosin staining after other staining.*

### 3. Other Staining

If immunofluorescence staining is carried out, after dyeing with eosin, wash the excess dye with distilled water and soak in tap water for 5min. Wash twice with 70% ethanol for 2min each time, then soak in PBS or physiological saline, TBS, TBST and other solutions for immunostaining or fluorescent dye staining for 5min, then immunofluorescent staining or other staining can be carried out.

## Note

1. The dye solution can be reused many times. When the effect is not good, can replace with the new dye solution. When there are a large number of samples, a dyeing frame and a dyeing cylinder can be used for easy operation.
2. When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test.
3. For your safety and health, please wear experimental clothes and disposable gloves.

## Related Products

G1010 Gimesa Stain Solution(Working Suit)

G1040 Wright Stain Solution

G1120 Hematoxylin-Eosin (HE) Stain Kit

G1140 Cole's Hematoxylin Solution (For Conventional Stain)

P2100 10×Polylysine