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规格: 2×50mL/2×100mL/2×250mL

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Cat: G1132

Size: 2×50mL/2×100mL/2×250mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	2×50mL	2×100mL	2×250mL	Storage
Reagent A: Fuchsin Solution	50mL	100mL	250mL	RT, avoid light
Reagent B: Methylene Blue Stain Solution	50mL	100mL	250mL	RT, avoid light

Introduction

The spores have high refraction, dense outer membrane, low permeability, and difficult to color and decolorize. When dyeing under heating condition, the dye not only enters into the cell but also into the spore. The dye entering into the cell can be decolorized by alcohol, while the dye in the spore still remains. After re-dyeing with counter stain, the cell and spore show different colors.

Self Provided Materials

95% Alcohol.

Protocol(for reference only)

1. Make a bacterial smear. Dry and heat-fix the smear .
2. Cover it with Reagent A. Heat the slide for 4-5 min; avoid boiling.If necessary, continue to drop the dye to avoid drying up. Maintained for 4-5 min.
3. Cool the smear. Decolorized with 95% alcohol until no red solution. Wash well with water.
1. Counterstain with Reagent B for 2-3 min. Wash with water until no dye comes out.
4. View under oil immersion lens.

Result

Thalli	Blue
Spore	Red

Note

1. The culture used for spore staining should be controlled for the age of bacteria, and most of spores should be kept in the body of bacteria.
2. When the dye is heated, it must be kept in the state of steaming out of the dye solution. The boiling will lead to the rupture of the cell or the spore sac. If the heating is not enough, the spore is difficult to be colored.
3. Please wear experimental clothes and disposable gloves.

Related Products

G1010 Giemsa Stain Solution (Working Suit)
G1020 Wright-Giemsa Stain Solution
G1060 Gram Stain Kit
G1160 Carbol-Fuchsin Solution