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Cole's Hematoxylin Solution (For Conventional Stain)

Cat:G1140

Size:100mL/500mL

Storage:RT, valid for at least 1 year.

Introduction

Hematoxylin is one of the most commonly used dyes in histochemistry and immunohistochemistry, which can be widely used for staining tissue sections or cultured cells. The staining solution can be used in combination with immunofluorescence staining or immunohistochemical staining. It can be stained with immunofluorescence or other dyes after hematoxylin staining solution, or it can be re-stained with hematoxylin after immunohistochemical staining. The dyeing process is simple and the operation time is short. The nucleus is blue after hematoxylin staining.

Protocol(for reference only)

Sample Treatment

For paraffin section: dewax in xylene twice for 5-10min each. Rehydrate with series ethanol (100%,90%,80%,70%) for each level 3min,and finally in distilled water for 2min.

For frozen section: restore to RT and wash with distilled water for 2min.

For cultured cell: fix with 4% PFA for more than 10min.Then wash with distilled water twice for 2min each.

Hematoxylin (HE) Staining

For the above treated samples, dye with Cole's Hematoxylin Solution (For Conventional Stain) for 5-10 min (the time can be adjusted according to the dyeing results and requirements), and then wash the remaining dye solution with tap water for about 10 min. Wash again with distilled water (for a few seconds).

Dehydration, Transparency, Sealing or Other Staining

a. Dehydration, transparency, sealing

- ① 95% ethanol twice each for 2 min.
- ② Transparent with xylene twice each for 5 min.
- ③ Seal with resinene or other sealing agent.
- ④ View under the microscope, the nucleus is blue.

b. Other Staining

If immunofluorescence staining is carried out, after dyeing with Hematoxylin, wash twice with 70% ethanol for 2min each time, then soak in PBS or physiological saline, TBS, TBST and other solutions for immunostaining or fluorescent dye staining for 5min, then can conduct immunofluorescent staining or other staining.

Note

1. When using this reagent for the first time, it is recommended to take 1-2 samples for pre-test.
2. When there are a large number of samples, a dyeing frame and a dyeing cylinder can be used for easy operation.
3. The dye solution can be reused many times. When the effect is not good, can replace with the new dye solution.
4. The staining process recommends light staining, usually only to distinguish the nucleus, too deep color may affect the color of the cytoplasm.
5. For your safety and health, please wear experimental clothes and disposable gloves.

Related Products

G1100 Eosin Y Stain Solution, For HE

G1120 Hematoxylin-Eosin(HE) Stain Kit

G1010 Gimesa Stain Solution(Working Solution)

G1040 Wright Stain Solution