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Cat: G1315

Size: 100mL/500mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

Neutral red acts as a pH indicator, changing from red to yellow between pH 6.8 and 8.0. The molecular formula is $C_{15}H_{16}N_4 \cdot HCl$, and the molecular weight is 288.78. In neutral or slightly alkaline environment, living cells of plants can absorb a large number of neutral red and excrete into vacuoles.

Neutral Red Solution, 0.5% can be used as a red nuclear counterstain in various histological procedures.

Protocol (for reference only)

1. Sample Treating

For Paraffin Section: Dewax in xylene twice for 10 min each. Absolute ethanol twice for 2 min each. Series ethanol (95%, 85%, 75%) for 2 min each level. Distilled water for 2 min.

For Frozen Section: Restore to room temperature for 2 min. Then wash with distilled water for 2 min.

For Culture Cell: Fix in 4% PFA for 15 min. Wash with distilled water twice and 5 min each time.

2. Neutral Red Staining

For the above treated samples, use Neutral Red Solution, 0.5% to stain for 0.5-1h.

Wash thoroughly with distilled water or tap water.

Viewing and photography (The dehydration process may cause discoloration, it is recommended to observe visually).

Note

1. If dehydration, transparency and sealing treatment are needed, xylene, Series Ethanol and resinene or other sealing agent shall be prepared by oneself.
2. When using this reagent for the first time, it is recommended to take 1-2 samples for pretest.
3. When staining with Neutral Red Solution, the time can be adjusted according to the dyeing results and requirements.
4. For your safety and health, please wear experimental clothes and disposable gloves.