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Mallory's Acid Hematoxylin (PTAH) Solution (Natural Oxidation)

Cat: G1381

Size: 100mL

Storage: RT, avoid light, valid for 2 years.

Introduction

Muscle fiber is a component of muscle tissue, which is composed of muscle cells. According to the morphological and functional characteristics, muscle fibers can be divided into smooth muscle (also known as striated muscle), skeletal muscle and cardiac muscle. There are many methods of muscle fiber staining, such as Ponceau method, Aniline Blue method, Hematoxylin Phosphotungstate method and so on.

Mallory's Acid Hematoxylin (PTAH) Solution (Natural Oxidation) is naturally oxidized and mature, with good staining ability and long preservation time. It is suitable for CNS, general tissue structure and all tissues fixed with standard fixative. Most of them are used to show the striae of striated muscle. This method is used to diagnose rhabdomyosarcoma. The staining time varies according to the preparation method, the fixative used and the displayed tissue structure. This product is only a single dyeing solution. Relevant auxiliary reagents should be prepared by yourself or refer to G1380-Mallory Acid hematoxylin dyeing kit (natural oxidation method).

Protocol(for reference only)

1. The tissue is recommended to be fixed in 10% formalin fixative or 4% PFA fixative, and conventional dehydration embedding.
2. Recommended paraffin section thickness 3-5 μ m, conventional dewaxing to distilled water.
3. (prepared by oneself) add the newly prepared PTAH oxidizing agent to cover the slice and oxidize it for 5min. Rinse with distilled water for 5-10s to remove excess reagents.
4. (Prepare by yourself) cover sections with oxalic acid solution and bleach for 1-2min. Rinse with tap water for 2min and balance with distilled water for 2min.
5. Immerse in Mallory's Acid Hematoxylin (PTAH) Solution (Natural Oxidation) and dye (cover) for 24h-48h.
6. Remove the slices and wash away excess liquid directly with 95% ethanol.
7. Dehydrate with anhydrous ethanol for 30 s, transparent with xylene for 1min, and seal with resinene.

Result

Cross-striations, fibrin, glial fibers, and nuclei	Deep Blue
Collagen Fiber, Cartilage Matrix	Brown Red
Coarse Elastic Fiber	Purple

Note

1. If the horizontal lines displayed after dyeing are not blue enough or bright red, it means that the dyeing liquid is insufficient or excessive oxidation, and the active ingredient content is insufficient, and it is necessary to change the dyeing liquid or prepare a new dyeing liquid.
2. This scheme avoids washing after dyeing, and also washes with 95% ethanol quickly to avoid background decolorization.
3. The dyeing range of this scheme will gradually increase with time, so overstaining should be avoided. After 24 h of dyeing, the staining degree can be observed under a microscope.
4. Dyeing can also be done at 60°C for 2-4 h instead of overnight dyeing at room temperature, but the effect may not be as good as room temperature dyeing.
5. For your safety and health, please wear a lab coat and disposable gloves.