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Nissl Stain Kit (Methyl Violet Method)

Cat:G1432

Size:2×50mL/2×100mL

Storage:RT, avoid light, valid for 6 months.

Kit Components

Reagent	2×50mL	2×100mL	Storage
Reagent(A):Methyl Violet Solution	50mL	100mL	RT, avoid light
Reagent(B): Nissl Differentiation	50mL	100mL	RT

Introduction

The neuronal soma consists of a large nucleus with a folded nuclear membrane, a sparse chromatin and a distinct nucleolus. In the cell body, the cytoplasm is Nissl granules, which can represent rough endoplasmic reticulum and produce specific punctate basophilic granules in many neurons. Nissl granules can be shown by many dyes, such as neutral red, methylene blue, toluidine blue and methyl violet. Variation of dyeing, pH and differentiation time make some staining not only highlight Nissl substance, but also include the nucleus and glia cells. Nissl body is a triangular or elliptical block substance distributed in the cytoplasm of nerve cells. It can be dyed purple-blue by methyl violet, methylene blue, toluidine blue and tar violet.

This kit has many characteristics: simple operation, stable staining, short differentiation time and wide application. It can be used in paraffin tissue sections of Nissl substance/Nissl bodies, neurons and so on.

Protocol(for reference only)

1. For paraffin section:fix fresh tissue in 10% neutral formalin fixative, then dehydrate and embed, the recommended thickness is 5-10 μ m; For frozen section:fix fresh tissue in 10% neutral formalin fixative, then dehydrate in series of sucrose for 24-72h and embed in OCT emding reagent, the recommended thickness is 10-15 μ m.
2. For paraffin section: dewax to distilled water. For frozen section: soak the section in distilled water and restore to room temperature.
3. Stain with Methyl Violet Solution for 10-20min.Rinse with distilled water to remove excess solution.
4. Differentiate with Nissl Differentiation for several seconds until the background is close to colorless.
5. Absolute alcohol dehydration, transparentize by xylene,seal with resinene.

Result

Nissl body	Violet Blue to Violet Black
Nucleus	Violet Blue
Background	Light Blue

Note

1. Nissl dissolves easily in vitro, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain. Tissue fixation plays a very important role. Ethanol, Carnoy or neutral formalin solution can be used for fixation.
2. This staining solution has a good effect on Nissl staining of paraffin tissue sections.
3. If only to confine the Nissl substance, it must be differentiated by acid differentiation solution after staining.
4. Paraffin section thickness is 5-10 μ m or 25 μ m (cortical neuron density should be assessed with 25 μ m section).
5. The stained specimens must be kept away from light, otherwise they will fade easily.
6. For your safety and health, please wear experimental clothes and disposable gloves.