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Nissl Stain Kit (Methylene Blue Method)

Cat: G1434

Size: 3×50mL

Storage: RT, avoid light, valid for 6 months.

Kit Components

Reagent	3×50mL	Storage
Reagent(A): Methylene Blue Stain	50mL	RT, avoid light
Reagent(B): Nissl Differentiation Solution	50mL	RT
Reagent(C): Ammonium Molybdate Solution	50mL	RT

Introduction

Nissl bodies are small triangular or oval shaped substances distributed in the cytoplasm of nerve cells, which can be dyed purple blue by basic dyes.

Nissl Stain Kit (Methylene Blue Method) has the advantages of simple operation, stable staining and wide application. It can be used for staining Nissl substances and neurons in paraffin sections. The presence and disappearance of Nissl body is an important indicator of whether nerve cells are damaged. When encephalitis, cerebral ischemia, axon reaction and other conditions occur, Nissl body will dissolve or even disappear.

Protocol(for reference only)

1. For paraffin section:fix fresh tissue in 10% neutral formalin fixative, then dehydrate and embed, the recommended thickness is 5-10 μ m; For frozen section:fix fresh tissue in 10% neutral formalin fixative, then dehydrate in series of sucrose for 24-72h and embed in OCT emding reagent, the recommended thickness is 10-15 μ m.
2. For paraffin section: dewax to distilled water. For frozen section: soak the section in distilled water and restore to room temperature.
3. Drop Methylene Blue Stain and stain for 10min.
4. Differentiate by Nissl Differentiation Solution for 1-3min(view under the microscope , stop differentiation until nissl body is clear).
5. (Optional)Treat in Ammonium Molybdate Solution for 3-5min.(See Note 3)
6. Rinse quickly in distilled water to avoiding decolorizing.
7. Dehydrate quickly in absolute ethanol, transparent by xylene and seal with resinene.

Result

Nissl body	Blue
------------	------

Note

1. Nissl body is easy to dissolve in vitro, so the tissue should be fixed immediately after removal, otherwise it is difficult to stain.
2. Tissue fixation plays an important role. Ethanol, Carnoy fixative or neutral formalin can be used for fixation.
3. Ammonium molybdate treatment is used to enhance electron density and increase contrast. If the decolorization is severe or the effect is poor after ammonium molybdate treatment, this step can be omitted.
4. This staining kit has a good effect on Nissl staining of paraffin sections.
5. For your safety and health, please wear experimental clothes and disposable gloves.