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Cat: G1452

Size: 100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Introduction

Calcium is abundant in human body. As a scaffold to support human body, bone plays an important role in secretion, transportation, muscle contraction, nerve conduction and so on. Alizarin Red S is an anthraquinone derivative, which is sodium alizarin sulfonate. It can chelate with calcium carbonate or calcium phosphate to form orange red complex. Alizarin Red S can get more reliable results for a small amount of sediment. It is often combined with fast green or Mayer hematoxylin staining solution to form orange red precipitate, which is suitable for dyeing a small amount of calcium salt tissue.

Protocol(for reference only)

For paraffin section

1. Fix the tissue in 10% neutral formalin or ethanol, then dehydrate and embed.
2. Dewax the section to 95% ethanol.
3. Place the section vertically and air dry thoroughly.
4. Add Alizarin Red S Solution, 1%, pH 4.2 onto the section and stain for 5-10min.
5. Quickly rinse in distilled water to prevent excessive washing from causing discoloration and removal of calcium salts.(See Note 1)
6. (optional)Re-dyeing with counterstain solution and wash with distilled water for three times.
7. Dehydrate with anhydrous ethanol, transparent with xylene and seal with neutral gum.

For cultured cell

1. Remove the medium from the plate and wash twice with PBS.
2. Fix in 10% neutral formalin or 4% paraformaldehyde for 10-15min.
3. Discard the fixative and wash with distilled water for 3 times.
4. After the water is completely absorbed, add Alizarin Red S Solution, 1%, pH 4.2 slowly and dye for 5-10min.
5. Quickly rinse in distilled water to prevent excessive washing from causing discoloration and removal of calcium salts.(See Note 1)
6. Add a proper amount of anhydrous ethanol into each hole to prevent drying in the hole. Observe and photograph under the microscope.

Result

Calcium Deposits	Orange Red
------------------	------------

Note

1. The dyeing time of Alizarin Red S Solution should be determined according to the content of calcium salt. It should be observed under the microscope. When the calcium salt is deep orange red, take out and wash. If the dyeing time is too long, it is easy to appear dispersion phenomenon, generally for 5min.If you want to preserve calcium salt staining to the greatest extent possible, you can use anhydrous ethanol to wash the slices and cells.
2. After dyeing with Alizarin Red S Solution, the calcium deposits are birefringent.
3. When re-dyeing with fast green, the background is green. When re-dyeing with Mayer's hematoxylin staining solution, the nucleus is blue.
4. This method is particularly useful in identifying and detecting small amounts of calcium, such as abnormal calcium in the kidney (excessive calcium in the urine).
5. For your safety and health, please wear experimental clothes and disposable gloves.