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Alizarin Red S Staining Extraction Solution (CPC, 10%)

Cat: G1453

Size: 100mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Alizarin Red S is an anthraquinone derivative, which is a sodium salt of Alizarin Sulfonate. It can chelate with calcium salts in calcium carbonate or calcium phosphate to form an orange red complex, suitable for tissue sectioning and calcium salt deposition staining of cultured cells. For the convenience of analyzing calcium salt content and concentration, Alizarin Red S staining can be completely eluted with CPC reagent. The eluted solution is placed in a 96 well plate, and the calcium salt content is semi quantitatively analyzed by measuring the OD value at 560nm.

This product is suitable for analysis of calcium salt content of various Alizarin Red staining kits.

Protocol(for reference only)

Alizarin Red S Staining Extraction Solution (CPC, 10%) may precipitate after being stored for a certain period of time due to its high concentration, which is a normal phenomenon. It can be used after being heated and dissolved in a 37 °C water bath or a 60 °C oven.

1. Perform Alizarin Red S staining according to customer experimental requirements.
2. For tissue sections, add 400ul Alizarin Red S Staining Extraction Solution (CPC, 10%) onto the section and incubate at room temperature for 15-30 minutes; for cultured cells, six well plate add 1ml Alizarin Red S Staining Extraction Solution (CPC, 10%) and incubate at room temperature for 15-30 minutes, until the Alizarin Red S staining is completely washed off.
3. Collect the eluted Alizarin Red S solution using a centrifuge tube and shake well. Place in a 96 well plate, add 100uL to each well, with a total of three repeated wells, and measure the OD value at a wavelength of 560nm on an enzyme-linked immunosorbent assay (ELISA) reader.

Note

1. Due to human operation errors and machine accuracy errors, the same batch of samples should be tested on the same day as much as possible.
2. For your safety and health, please wear laboratory clothes and disposable gloves when operating.

Related products

G1450 Alizarin Red S Solution, 0.2%, pH8.3

G1452 Alizarin Red S Solution, 1%, pH4.2

G3280 Calcium Stain Kit (Modified Alizarin Red S Method)

G3281 Calcium Stain Kit (Alizarin Red S Method)

G3283 Alizarin Red S Staining Quantitative Detection Kit (CPC Method)