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## Elastic Fiber Stain Kit (Orcein Method)

**Cat:** G1589

**Size:** 3×50mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

### Kit Components

Reagent	3×50mL	Storage
Reagent(A): Orcein Solution	50mL	2-8°C
Reagent(B): Mayer Hematoxylin Solution	50mL	RT, avoid light
Reagent(C): Scott bluing Solution	50mL	RT

### Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elasticize usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin-fuchsin, and Weigert's elastic stain in histological sections.

Elastic Fiber Stain Kit (Orcein Method) is a selective staining method for identifying elastin protein. Sometimes the lipofuscin in the hepatocytes will be stained.

### Protocol(for reference only)

1. Fix the tissue in 10% neutral formalin fixative, and routinely dehydrate and embed. Cut the slice in 4-5μm and dewax to 70% alcohol.
2. Soak in Orcein solution and stain for 6-8h at 37°C incubator or stain for 12-24h at room temperature.
3. Rinse in 70% alcohol for 3-5s, and then wash with distilled water for 1-2min.
4. (Optional)Re-dyeing with Mayer Hematoxylin Solution for 2-3min.
5. (Optional)Rinse in distilled water for 2-3min, differentiating in 1% acid alcohol if necessary.
6. (Optional)Rinse in Scott blue Solution for 2min, and wash with distilled water.
7. Dehydrate with series alcohol, transparent with xylene and seal.

### Result

Elastic Fiber	Brownish Red to Dark Brown
Background	Grayish Brown or Reddish in Varying Degrees

### Note

1. The Orcein Solution should be restored to room temperature before use, and it can be reused several times.
2. The dyeing time of elastic fibers can be appropriately shortened or extended according to different tissues.
3. For your safety and health, please wear experimental clothes and disposable gloves.