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Cat: G1594

Size: 100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Introduction

In a narrow sense, connective tissue contains three types of fibers: collagen fibers, reticular fibers, elastic fibers. Elastic fibers are found in the lungs, arteries, veins. It is highly refractive, elasticized usually thinner than collagen fibers. Elastic fibers stain well with Gomori aldehyde fuchsin, orcein, resorcin- fuchsin, and Weigert's elastic stain in histological sections.

Aldehyde-Fuchsin Stain Solution can be used for many purposes. Depending on the staining time, it can be used for staining elastic fibers, as well as for staining lipofuscin, mast cells, and pancreatic islets β cells and pituitary basophils.

Protocol (for reference only)

1. Fix in 10% neutral formalin fixative, and routinely dehydrate and embed.
2. Cut the paraffin embedded tissue into 4 μ m-7 μ m thin sections and routine dewax to water.
3. Operate according to the specific experimental requirements. If dyeing elastic fibers, please refer to the steps in the G1593 manual.

Result

Elastic Fiber	Purple to dark purple
Mast cell granules and mucus	Purple to dark purple
Background	Different degrees of yellow

Note

1. This method can display elastic fiber, anterior elastic fiber and acid resistant fiber, but it needs to slice slightly thicker, 7 μ m is the best.
2. When aldehyde fuchsin is dyed, it should be capped to prevent the solution from volatilizing.
3. When the Aldehyde-Fuchsin Stain Solution is stored for a long time, the dyeing power will decrease, and the dyeing time should be increased.
4. When islet beta-cells were stained, the time should controlled at 30 min; when pituitary basophils were stained, the time is controlled at 60 min.
5. For your safety and health, please wear experimental clothes and disposable gloves.