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规格: 4×100mL/4×500mL

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**Cat:** G1617

**Size:** 4×100mL/4×500mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	4×100mL	4×500mL	Storage
Reagent(A):Nuclear Staining Solution	100mL	500mL	RT, avoid light
Reagent(B):Cytoplasm Staining Solution	100mL	500mL	RT, avoid light
Reagent(C):Decolorization Solution	2×100mL	2×500mL	RT

### Introduction

The Papanicolaou method is commonly used for routine cytological staining. Papanicolaou stain was initially only used to detect vaginal epithelial estrogen levels and pathogens such as Candida and Trichomonas in the genital tract. Pasteurized staining solution contains a variety of ions, which can combine with proteins with opposite charges in cells, so as to distinguish and display various components in cells. After staining, the cytoplasm is bright, transparent, and the nuclear structure is clear. The cytoplasm of keratinized squamous cells is pink or orange, the non keratinized cells are green, and the nucleus is blue or black.

This product is a non-toxic rapid Papanicolaou staining kit carefully improved by our company. Compared with the classical PAP staining method, the steps are cumbersome and time-consuming, which simplifies the operation steps and dyeing time. The whole dyeing process is only divided into three steps, which can be completed in 5min, without gradient alcohol dehydration and xylene transparency, and can be directly sealed and stored. The kit is mainly used for exfoliative cell staining, which is applicable to various specimens of epithelial cells and mesothelial tissues. It is the most commonly used staining method in vaginal exfoliative cell examination. It can also be used for cervical scraping, leucorrhea, pleural effusion and ascites, sperm, TCT and other samples.

### Self Provided Materials

AF Fixative(alcohol-formalin fixative), series of ethanol, optical microscope, resinene.

### Protocol(for reference only)

1. Take 5-10ul of the sample to be examined and evenly coat it on the anti peeling section to prepare the cell smear.
2. After the cell smears naturally dry, fix with AF fixative or 95% ethanol for 2-3min and wash with distilled water for 30s.
3. Drop stain with Reagent (A): Nuclear Staining Solution for 1-2min, rinse with tap water for 30-60s and return to blue.
4. Drop stain with Reagent (B): Cytoplasm Staining Solution for 1-2min.
5. Rinse with Reagent (C): Decolorization Solution for 5-10s until no excess dye occurs, then dry the smear naturally.
6. (Optional)Dehydrate by series of ethanol and transparent by xylene.
7. Seal with resinene.

### Result

Nucleus	Blue Purple or Black
The cytoplasm of non keratinized cells	Light Blue or Light Green
The cytoplasm of keratinocytes	Pink Red or Orange Red

### Note

1. The cleaning steps throughout the entire process need to be handled gently to prevent cell shedding.
2. The staining solution that stored for a long time may cause sedimentation, which is a normal phenomenon and can be filtered before use.
3. For your safety and health, please wear laboratory clothes and disposable gloves when operating.