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Fibrin Stain Kit (Modified MSB Method)

Cat: G2040

Size: 9×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	9×50mL	Storage
Reagent (A): Hypo Solution	50mL	RT
Reagent (B): Celestite Blue Solution	50mL	2-8°C, avoid light
Reagent (C): Mayer Hematoxylin Solution	50mL	2-8°C, avoid light
Reagent (D): Acid Differentiation Solution	50mL	RT
Reagent (E): Mathew Yellow Solution	50mL	RT, avoid light
Reagent (F): Acid Red Solution	50mL	RT, avoid light
Reagent (G): Phosphotungstic Acid Solution	50mL	RT, avoid light
Reagent (H): Aniline Blue Solution	50mL	RT, avoid light
Reagent (I): Weak Acid Solution	50mL	RT

Introduction

The pathological endogenous precipitate is a part of pigment precipitate. After certain pathological changes, the tissue cells form the precipitate in different shapes and characteristics. This special protein polymer can show cellulose protein after dyeing. Cellulose is a special protein formed by the aggregation of fibrin molecules in the blood, also known as fibrin. Under normal circumstances, this kind of protein exists in the tissue in the form of curved filamentous cellulose, most of which are in the form of network structure, sometimes in the form of thick cellulose network, and old agglutinable in the form of amorphous block. When the tissue is damaged, the vascular endothelium is seriously damaged, and the vascular permeability increases accordingly, which can lead to the leakage of a large number of fibrin.

This kit improved based on the classic MSB method, which takes Mathew Yellow Solution, Acid Red Solution and Aniline Blue Solution as the core, which is composed of Celestite Blue Solution, Mayer Hematoxylin Solution, Mathew Yellow Solution and Aniline Blue Solution. It is a simple and cheap dyeing solution. After dyeing, the color of cellulose is red or blue.

Protocol (for reference only)

1. Conventionally cut paraffin section and conventionally dewax.
2. Treat in Hypo Solution for 3-5mins, wash with distilled water to remove excess solution.
3. Stain in Celestite Blue Solution for 3-5mins. Wash with distilled water for 1min.
4. Stain in Mayer Hematoxylin Solution for 3-5mins. Wash with distilled water for 1min.
5. Differentiate by Acid Differentiation Solution for several seconds, then rinse with tap water for 10mins.
6. Stain in Mathew Yellow Solution for 2-3mins and quickly wash with distilled water for 5-10s.
7. Stain in Acid Red Solution for 10mins and quickly wash with distilled water for 5-10s.
8. Treat the section with Phosphotungstic Acid Solution until the collagen red disappears, and quickly wash with distilled water for 3-5s.
9. Stain in Aniline Blue Solution for 3-5mins or more. Slightly wash with Weak Acid Solution.
10. Quickly dehydrate by 95% ethanol for 3-5s, absolute ethanol twice and each for 5-10s.
11. Transparent by xylene twice and each for 1-2mins. Seal with resinene.

Result

Cellulose	Red (Fresh cellulose can be yellow, old cellulose can be blue)
Nucleus	Blue-Gray
Erythrocyte	Yellow
Collagen Fiber	Blue

Note

1. Phosphotungstic Acid Solution treatment generally takes 5-10mins, at least more than 5mins to produce.
2. Aniline Blue Solution staining is usually for 5mins, but it should be observed every 2mins to prevent over staining, otherwise the blue is too deep.
3. For your safety and health, please wear experimental clothes and disposable gloves.