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Cat: G3060

Size: 100mL

Storage: -20°C, avoid light, valid for 6 months.

Kit Components

Reagent	Size	Storage
Reagent(A): GUS Buffer A	50mL	RT
Reagent(B): GUS Buffer B	0.2mL	2-8°C, avoid light
Reagent(C): GUS Buffer C	0.2mL	2-8°C, avoid light
Reagent(D): 2×Fixation Buffer	25mL	RT, avoid light
Reagent(E): X-GlcA	80mg	-20°C, avoid light
Reagent(F): X-GlcA Solvent	1mL	RT, avoid light
Reagent(G):Deionized Water	100mL	RT

Introduction

X-Gluc(X-GlcA) has a molecular weight of 521.8 and CAS number of 114162-64-0. It is a substrate for the detection of GUS gene in *Bacillus coli* and can rapidly detect GUS gene fusion markers in plants. Under suitable reaction conditions, β -glucosidase (GUS) can hydrolyze X-Gluc into blue substance, which is not dissolved in Indigo substance in transgenic nuclear tissue. The sites or sites with GUS activity show blue or blue spots, which can be observed by naked eye or microscope. The GUS Staining Kit is mainly used for GUS gene expression analysis of transgenic plants.

Protocol(for reference only)

1. Prepare X- GlcA Solution: add 229ul of X-GlcA Solvent to 80mg X-GlcA, or take appropriate amount of the above two substances to dissolve, so that the concentration reaches 350mg/mL, mix them gently to form X-GlcA Solution. After distributing, store them in dark at - 20 °C.
2. Prepare GUS Solution according to the radio as follows.

Reagent	volume
GUS Buffer A	2.5mL
GUS Buffer B	10ul
GUS Buffer C	10ul
Deionized water	5.5mL
Methanol	2.0mL
X-GlcA Solution	20ul
Total volume	Approximately 10mL

3. Fixation(Fixation is not a necessary step. If fixation is not required, proceed to step 4):
 - ①Take the transgenic plant tissue and add it into a about 3-5mL of clean vial or porous plate.
 - ②Take equal amount of 2×Fixation Buffer and Deionized Water to mix to form 1×Fixation Buffer. Add 1×Fixation Buffer to completely submerge the tissue. Incubate at room temperature for 45min, discard the liquid.
 - ③Prepare Washing Solution: mix GUS Buffer A with Deionized Water according to the radio of 1:20 to form Washing Solution. Rinse 3 times with the prepared Washing Solution for 1 min each time.
4. Add appropriate GUS Solution to completely submerge the tissue.
5. Use a vacuum pump to extract the gas in a vial or porous plate, and the extraction time shall be more than 2min. The purpose of this step is to extract gas from plant tissues and make the dye more easily enter the tissues.
6. Immediately cover the vial or porous plate and incubate at 37 °C for 24h. With the increase of incubation time, blue gradually appears. When the expression level is high, the sites or sites of GUS activity show blue or blue spots.

7. Remove the chlorophyll of the sample with 75% ethanol, and generally soak the sample in ethanol for 1-3h. If necessary, repeat the decolorization step to completely remove chlorophyll. The samples are stored in ethanol and can be observed by naked eye or ordinary optical microscope.

Note

1. The prepared GUS Solution can be stored for half a month at 2-8 °C in dark.
2. X-GlcA Solution should avoid repeated freezing and thawing, otherwise the dyeing efficiency will decrease.
3. Due to the tissue specificity and other reasons, the blue color reaction may not be completely consistent, so pay attention to explore the specific experimental conditions.
4. For your safety and health, please wear experimental clothes and disposable gloves.