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**Cat:** G3225

**Size:** 5×50mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit Components

Reagent	5×50mL	Storage
Reagent A: Crystal Violet Staining Solution	50mL	RT, avoid light
Reagent B: Lugo iodine solution	50mL	RT, avoid light
Reagent C: Compound red staining solution	50mL	RT, avoid light
Reagent D: Differentiation fluid	50mL	RT
Reagent E: Counterstain	50mL	RT, avoid light

### Introduction

The Gram stain, invented in 1884 by the Danish physician Christain Gram, is a widely used identification stain in bacteriology. By staining with this method, bacteria can be identified into two major groups, Gram-positive (G+) and Gram-negative (G-). The different coloring reactions of bacteria are due to differences in cell wall permeability to ethanol and resistance to decolorization, mainly determined by the thickness and structure of the peptidoglycan layer.

This staining kit is modified from the original Gram staining method and is particularly suitable for identifying bacterial infections in tissue sections. Bacteria are stained dark purple or pink, and the background of the tissue is colored yellow in the presence of a re-staining solution, allowing a clear view of the bacteria present in the tissue.

### Protocols(for reference only)

1. Use paraformaldehyde to fix the sections normally, the embedded tissues were cut into 3-7 $\mu$ m thin slices, deparaffinized to distilled water, and then the following staining steps were performed. Crystalline violet staining solution was used to stain the sections for 2 min, rinsed with tap water, and excess staining solution was removed.
2. Crystalline violet staining solution was stained for 2min, rinsed with tap water to remove excess staining solution.
3. Mordant staining with Lugol's iodine solution for 5min, rinse with tap water, remove excess iodine, drain off excess water but not completely dry.
4. Use 95% ethanol to differentiate until no excess purple elution, rinse with tap water to terminate differentiation.
5. Stain with red staining solution for 5 min, wash briefly with tap water to remove excess staining solution.
6. Add the differentiation solution dropwise for 5min, rinse with tap water, drain off excess water but not completely dry.
7. 95% ethanol rapid wash three times dehydration.
8. Rapid staining with re-staining solution for 10-20s, washing with anhydrous ethanol for 5-10s.
9. Anhydrous ethanol-xylene (anhydrous ethanol and xylene 1:1 mixture) dehydration 5-10s.
10. Xylene transparent, neutral gum sealing.

### Result

Gram-positive bacteria and nuclei	Purple
Gram-negative bacteria	Red
Background	Yellow

### Note

1. The bacteria are small and the results need to be observed with a 100x oil microscope.
2. 95% ethanol differentiation step should be replaced with new 95% ethanol several times until no purple color comes off if the drop is differentiated.
3. For your safety and health, please wear lab coat and disposable gloves to operate.