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Luxol Fast Blue Myelin Stain Kit(Cresyl Violet Method)

Cat: G3245

Size: 3×50mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent	3×50mL	Storage
Reagent(A): LFB Staining Solution	50mL	RT, avoid light
Reagent(B): LFB Differentiation Solution	50mL	RT
Reagent(C): Cresyl Violet Staining Solution	50mL	RT, avoid light

Introduction

Myelin sheath (myelin sheath) is a myelinated nerve fiber axon externally wrapped tubular outer membrane, mainly composed of myelin, myelin sheath has Langfeld's junction, which allows nerve impulses to jump conduction. Myelin staining has certain significance in pathological diagnosis. The pathological changes of myelin sheath can be divided into early, middle and late stages. In the early stage of the lesion, the color is deep. In the middle stage, the myelin degenerated into lipid droplets. In the late stage, the myelin sheath degenerated completely and is cleared by phagocytes, so there is no positive result.

Luxol Fast Blue Myelin Stain Kit can show whether the myelin sheath is complete, denatured, necrotic and repaired under pathological conditions. It has significance for the pathological diagnosis and research of nerve tissue.

Protocols(for reference only)

1. Cut into paraffin section in 5-8 μ m thick, then dewax to 95% ethanol.
2. Sections were immersed into LFB Staining Solution and stained overnight at room temperature (paraffin sections are recommended for 12-18h, frozen sections were stained for no more than 16h).(See Note 4)
3. Wash the excess dye solution in 95% ethanol, rinse in distilled water.
4. Decolorize with LFB Differentiation Solution for 3-5s, and pour off the differentiation solution.
5. (optional)Decolorize with 70% ethanol for 30s until the gray and white matter was clear.
6. Rinse in distilled water(if the color differentiation is insufficient, can repeat steps 4 and 5).
7. Re-dyeing with Cresyl Violet Staining Solution for 30-40s, wash with distilled water.
8. Immediately dehydrated rapidly with 95% and 100% ethanol, transparent by xylene clear and seal with resinene.

Result

Myelin Sheath	Blue
Neuron	Pink Red to Purple

Note

1. Differentiation is a key step. The differentiation time should be strictly controlled, and the degree of differentiation can be observed under the microscope.
2. It is better to use 10% formalin as fixative.
3. The slices should not be too thick, paraffin sections should be controlled within 3-7 μ m, frozen sections should be controlled within 10-20 μ m, otherwise it is easy to appear phenomena such as dehiscence or overstaining.
4. If the dyeing effect at room temperature is not good, can stain at 56 °C for 2h.
5. If the color separation of the LFB Differentiation Solution is too strong, it can be diluted 2-3 times with distilled water before use, while paying attention to the differentiation time.
6. For your safety and health, please wear experimental clothes and disposable gloves.