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Modified Nerve Silver Plating Staining Kit, Bielschowsky's Method

Cat: G3260

Size: 5×50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

Reagent	5×50mL	Storage
Reagent(A): Pretreatment Solution	50mL	2-8°C, avoid light
Reagent(B): Chromogenic Solution	50mL	RT
Reagent(C): Ammoniac Silver Solution	50mL	2-8°C, avoid light
Reagent(D): Toner Solution	50mL	2-8°C, avoid light
Reagent(E): Hypo Solution	50mL	RT

Introduction

Neurons, also known as nerve cells, are the basic units that make up the structure and function of the nervous system. A neuron consists of a cell body and a protrusion, which can be further divided into axons and dendrites. Myelinated nerve fibers have a sheath over the axon called myelin, and the tiny branches at the end of it are called nerve endings, which can form motor endplates with target muscles. There are many staining methods for neurons and nerve fibers, such as silver staining and cresyl violet staining.

Modified Bielschowsky's Stain Kit is a typical silver staining method. Its basic principle is that the fixed tissue and sections are immersed in the silver solution, and then treated with reductant, so that the silver particles settle in the axonal plasma and appear dark brown or black. After silver staining, many interlaced filaments can be seen in the cytoplasm of neurons, and extend to dendrites and axons. The staining method shows that neurofibroma and ganglioneurofibroma are positive, while schwannoma is negative.

Protocols (for reference only)

The Chromogenic solution is stock solution. Before use, dilute 4 times with distilled water to form Chromogenic Working Solution. (See Note 3)

1. Cut paraffin section in 8-15 μ m thick, dewax to distilled water.
2. Add the section in Pretreatment Solution, and incubate in a 37°C incubator for 25-35min in dark. Rinse in distilled water for 2-3min.
3. Reduce with Chromogenic Working Solution for 10-20s until the section turns yellow. (See Note 4)
4. Rinse in distilled water for 3-5min.
5. Drop Ammoniac Silver Solution onto the section and stain for 1-2min.
6. Add equal amount of Chromogenic Working Solution on the slice, and incubate in dark for 2-10 min.
7. After pouring, add the Chromogenic Working Solution again to reduce for 1 min and fully react. (See Note 5)
8. Treat with Hypo Solution for 3-5min. Rinse in distilled water for 1-2min.
9. Match color with Toner Solution for 1-5min.
10. Rinse in distilled water for 3-5min, then use filter paper to absorb the moisture around the section.
11. Dehydrate in 95% ethanol and absolute ethanol, transparent by xylene and seal with resinene.

Result

Neurons, Axons, Nerve Fibers	Dark Purple to Black
------------------------------	----------------------

Note

1. The glass container used should be very clean. Wash repeatedly with water and distilled water.
2. In silver staining, the sections should be flattened to avoid wrinkles, so as to avoid uneven coloring.
3. The Chromogenic Solution is a reserve solution, and direct use will cause non specific and uneven background adhesion. It is recommended to dilute it 4 times with distilled water to make a working solution for use.
4. Frozen sections can be observed to turn yellow, while paraffin sections may not have significant discoloration reactions, which is a normal phenomenon.
5. Directly adding a reducing agent to the Ammoniac Silver Solution can to some extent avoid non-specific coloring. Under normal circumstances, the nerve fibers are copper red to black in color, with no obvious particle deposition. After dyeing the slices, it is important to handle them gently when mounting them to avoid crushing them.
6. For your safety and health, please wear laboratory clothes and disposable gloves for operation.