## Dacie 氏液

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**Storage:** RT, valid for 6 months.

### Introduction

Red blood cell (RBC) is the most numerous blood cell in the blood. The principle of Dacie's fluid is to use isotonic diluent to dilute the blood in a certain number of times, fill it into the counting chamber, count the number of RBC in a certain volume under the microscope, and calculate the number of RBC in each liter of blood by conversion. The Dacie's Solution is only used in scientific research field, not in clinical diagnosis.

### Self Provided Materials

Fresh whole blood, Micropipette, Cell counting plate, Microscope

### Protocol(for reference only)

1. Take small test tube and add 4.0mL Dacie's Solution.
2. Take 20 $\mu$ L of peripheral blood with a clean and dry micropipette, wipe off the remaining blood outside the micropipette, add it to the bottom of Dacie's Solution, gently drain the blood, then gently suck the upper supernatant to clean the pipette for 2-3 times, and mix it immediately.
3. After mixing well, fill the red blood cells with a clean micropipette into the counting chamber. Pay attention to prevent bubbles or overflows. Leave it at room temperature for 1-2 min.
4. Under the high power microscope, count the platelets in the four corners of the central square grid and the central 5 square grids. Pressure line cells are counted according to the principle of "counting up and not counting down, counting left and not counting right".

### Calculation

The number of red blood cell/L = The number of red blood cells in 5 middle squares $\times 5 \times 10 \times 200 \times 10^6$   
=The number of red blood cells in 5 middle squares $\times 10^{10}$   
=The number of red blood cells in 5 middle squares/ $100 \times 10^{12}$

$\times 5$	5 middle squares convert to 1 large square
$\times 10$	The volume of a square grid is 0.1 $\mu$ L, which is converted into 1.0 $\mu$ L.
$\times 200$	The actual dilution ratio of blood shall be 201 times, but 200 times for calculation.
$\times 10^6$	Convert 1 $\mu$ l to 1L

### Note

1. Avoid squeezing too much when collecting blood, and the acupuncture depth should be appropriate.
2. The small test tube and counting plate should be clean to avoid misidentification of cells.
3. Within the reference range, the difference between the two counts of red blood cells shall not exceed 5%.
4. Hemoglobin concentration is not allowed to convert the number of red blood cells.
5. For your safety and health, please wear experimental clothes and disposable gloves.