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Mast Cells Stain Kit(Alcian Blue-Safranin Method)

Cat: G3672

Size: 3×50mL

Storage:RT, avoid light, valid for 6 months.

Kit Components

| Reagent | 3×50mL | Storage |
|---|--------|-----------------|
| Reagent(A): Csaba Acidizing Solution | 50mL | RT |
| Reagent(B): Alcian Blue-Safranin Solution | 50mL | RT, avoid light |
| Reagent(C): Csaba Dehydrating Solution | 100mL | RT, avoid light |

Introduction

Mast cells are common cells in loose connective tissue. They are often distributed in groups along small blood vessels and lymphatics. They are also common around interlobular ducts of bronchi and pancreas. Generally, mast cells are large, with a diameter of about 20-30 μm , round or oval. Their nucleus are small, and their cytoplasm is filled with thick and heterochromatic eosinophils.

Mast Cells Stain Kit(Alcian Blue-Safranin Method) can distinguish different types of metachromatic particles. After staining, immature mast cell particles are blue, mature mast cell particles are red, which is relatively clear.

Self Provided Materials

10%NBF, Series Ethanol, Dye Vat

Protocol(for reference only)

1. Fix the fresh tissue immediately in 10% NBF,conventional dehydration and embedding.
2. Cut into 5-6 μm thick sections, and dewaxing to water.
3. (Optional)Add Csaba Acidizing Solution for 3mins.
4. Stain with Alcian Blue-Safranin Solution for 15-30mins,Rinse with distilled water.
5. Dehydrate with Csaba Dehydrating Solution repeated.(Control under microscope)
6. Transparent with xylene and seal with resinene.

Result

| | |
|--------------------|------|
| Immature Mast Cell | Blue |
| Mature Mast Cell | Red |

Note

1. For mast cells staining, the tissue should be fresh and fixed immediately after removal.
2. If there is no 10% NBF, it can be replaced by normal saline containing 4% FA.
3. If the Csaba Dehydrating Solution is empty,ethanol should be used as an alternative. By the way,the section may fade easily.
4. After differentiation, the time of transparent should be not for a long time, otherwise it is easy to cause fading.
5. For your safety and health, please wear experimental clothes and disposable gloves.