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Plant Tissue ROS Detection Kit(DAB Method)

Cat: G4815

Size: 3×100mL

Storage: 2-8°C, avoid light, valid for 1 year.

Kit Components

	Reagent	3×100mL	Storage
Reagent A: DAB Working Solution	Reagent A1:DAB Store Solution	10mL	2-8°C, avoid light
	Reagent A2:DAB Diluent Solution	90mL	2-8°C
Before use, mix A1 with A2 by 1: 9 to prepare DAB Working Solution.			
Reagent B:Decolorization Solution		100mL	RT
Reagent C:Preservation Solution		100mL	RT

Introduction

Plant tissues produce various reactive oxygen species (ROS) under stress environmental conditions, and ROS activity is very high and extremely unstable. Therefore, the detection of ROS is usually determined by its final product. Peroxides are one of the main components of plant reactive oxygen species. Under the catalysis of peroxidase, peroxides can quickly react with DAB (3,3-diaminobenzidine tetrahydrochloride) to form brownish red precipitates, thereby locating peroxides in tissues.

The Solarbio Plant Tissue ROS Detection Kit(DAB Method) is used for peroxide staining in plant living tissues. Generally used for overall staining of tender root tips, leaves, etc. After staining, the areas in the tissue where peroxides accumulate appear brown to dark brown.

Protocol(for reference only)

- Organizational Preparation:**Collect plant seedlings or leaves that have been treated (such as heavy metal stress), wash them with tap water, and place them on filter paper to absorb excess water.
- Tissue Staining:**Immerse the experimental sample in DAB Working Solution, treat it under negative pressure of -0.1MPa for 30 min, and then immerse it in room temperature and dark light for 4-12 hours until the positive area appears dark blue, and the remaining areas appear light blue or nearly colorless or in the color of the plant itself. (Negative pressure and staining time can be adjusted based on the degree of plant tenderness and color development)
- Tissue Decolorization:**Carefully remove the experiment with tweezers, immerse it in distilled water and rinse back and forth 3-5 times. Place it on a filter paper to absorb excess water, immerse it in the tissue decolorization solution, and treat it in a water bath at 70-80 °C for 20-40 min until the background color of the tissue is completely removed. If the color of the decolorization solution is darker during the treatment, replace it with fresh tissue decolorization solution.
- Result Observation:**After cooling, take out the experimental sample, immerse it in distilled water and rinse it back and forth 3-5 times. Place it on a filter paper to absorb excess water, then transfer the sample to an appropriate amount of tissue preservation solution and soak it for 10 to 30 min. Then, take out the sample and take a photo. The sample can be placed in a storage solution and stored at room temperature for more than a month.

Result

ROS Positive Areas	Brown to Red Brown
Background(after decolorization)	Colorless or Light Yellow

Note

- Peroxides are prone to decomposition, so plant samples need to be collected fresh and stained as soon as possible.
- Any external stimuli may stimulate plant stress to produce peroxides, and it is recommended to take complete samples as much as possible to avoid artificial damage causing false positives.
- It is usually recommended that the volume of the decolorization solution be more than 10 times the sample volume during decolorization. It can be reused but it is recommended to discard it when the color is too dark.
- After the staining of the tissue sample is completed, it is necessary to take photos and save the results as soon as possible.
- For your safety and health, please wear laboratory clothes and disposable gloves for operation.