

Helicobacter pylori selective additive

Cat: LA7471

Specification: 10*1ml

Storage: Store at 2-8°C, and it is valid for 1 year.

Product Information:

The dose of vancomycin	0.001g
Amphotericin B	0.5mg
Trimethoprim	2mg
Cefsulodin	1mg
distilled water	1ml
property	colorless transparent

Instructions for use: (for reference only)

In a sterile environment, use a grinding wheel to cut the neck of the ampoule bottle directly below the blue point, break the ampoule bottle away from the blue point (wear gloves), when the medium is cold to 45-50°C, use a pipette to add appropriate amount of this product into the medium during aseptic operation, mix gently and set aside.

Note:

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.

Related Products:

LA7470 *Helicobacter pylori* AGAR base

